Photometer PF-12

Handbuch
Manual
Manuel

Kompaktphotometer PF-12
Compact photometer PF-12
Photomètre compact PF-12
Fotómetro compacto PF-12

MACHERY-NAGEL

www.mn-net.com

Deutsch	3
English	17
Français	31
Español	45

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CONTENTS

1. Safety precautions	18
1.1. General safety precautions	18
1.2. Chemical and biological safety	18
2. Photometer PF-12	19
2.1. Technical Data	19
2.2. Design	20
2.3. Shipment	20
2.4. Interface and power supply	20
3. Usage	21
4. Setup	22
4.1. Unpacking the instrument	22
4.2. Operation with batteries	22
4.3. Operation with battery pack	22
5. Operation	23
5.1. Keyboard	23
5.2. Configuration	23
5.3. Performance of <i>VISOCOLOR® ECO</i> tests	24
5.4. Performance of <i>NANOCOLOR®</i> tube test	25
5.5. Basic photometric functions	26
6. Memory	27
7. Special functions	28
7.1. Special methods	28
7.2. Sample numbering	29
7.3. Identification of sample location	29
7.4. Sample dilution	30
8. Cleaning the photometer	30
9. Trouble shooting	30

1. Safety precautions

1.1. General safety precautions

Read the entire manual carefully before you set-up and start the instrument. Please adhere to all information. Non-compliance may result in injury to the users or damage to the instrument. In order to not adversely influence the instrument's inbuilt safety measures, the instrument may not be installed or used in any other way than specified in the in the manual.

1.2. Chemical and biological safety

Potential dangers can arise due to exposition to chemical or biological substances.

Working with chemical samples, standards and reagents bears certain risks. We recommend for all users of such products to look into safety measures, proper use of chemicals and the respective MSDS prior to using the instrument and chemicals.

When using the instrument normally, it may be necessary to handle dangerous substances or biologically harmful samples.

- Before handling any such dangerous / harmful substances, you need to pay attention to all hazard statements and safety precautions on the chemical's package and in the MSDS.
- All used solutions have to be disposed of according to national law and regulations.
- Chose your safety equipment to fit the concentration and amount of hazardous chemical at your work place.

2. Photometer PF-12

2.1. Technical Data

Type:	Single-beam filter photometer with microprocessor control, self-test and auto-calibration
Wavelength range:	340–860 nm
Optics:	Automatic filter wheel with 7 interference filters, Insensitive to external light for fast measurements without cuvette slot cover
Wavelengths:	345 / 436 / 470 / 540 / 585 / 620 / 690 nm plus 1 compartment for an additional filter
Wavelength accuracy:	± 2 nm Bandwidth at half of maximum transmission 10–12 nm
Light source:	Tungsten lamp
Detector:	Silicon-photodiode
Blank value:	Automatic
Measuring modes:	Over 100 preprogrammed tests (<i>NANOCOLOR</i> ® tube tests und <i>VISOCOLOR</i> ® <i>ECO</i> tests), absorbance, transmission, factor, standard, 20 freely programmable methods
Photometric range:	± 3 E
Photometric accuracy:	± 1 %
Stability:	< 0.002 E/h
Cuvette holder:	Round tubes 16 mm OD
Data memory:	200 results, GLP-compatible
Display:	Backlit graphic display, 64 x 128 pixels, 12 languages; All important data at a glance: result with measuring unit, date, time, sample number, sample location, dilution
Operation:	Intuitive menu navigation, durable keypad, test selection via entering the test number or from parameter lists
Quality control:	With <i>NANOCONTROL</i> NANOCHECK
Interfaces:	USB 2.0
Update:	Free of charge via Internet / PC
Operating range:	0–50 °C, up to 90 % relative humidity
Power supply:	USB-power supply, (rechargeable) batteries, battery pack
Housing:	Waterproof, IP 67
Dimensions:	215 x 100 x 65 mm
Weight:	0.7 kg
Warranty:	2 years

This device complies with the following directives:

2006/95/EC Low Voltage Directive

2004/108/EC EMC Directive

2.2. Design

1. Graphic display 128 x 64 Pixels
2. 20-part keypad
3. Cuvette slot for tubes with 16 mm outer diameter
4. USB-interface (2.0) for connection of power supply, charging the optional battery pack (as of SN: NPF121600), data transfer and update

2.3. Shipment

If one of the components from the following list should be defective or missing, please contact MACHEREY-NAGEL or your local distributor.

Compact photometer PF-12 incl.:

Software DVD, manual, 4 AA Mignon batteries, 4 empty cuvettes, funnel, beaker, syringe, USB-cable, calibration cuvette and certificate in in a rugged case

2.4. Interface and power supply

USB-interface for:

- Update via Internet
- Data export to PC
- Power supply via PC USB-interface
- External power supply (REF 919 220) with USB-interface

Additional power supply:

- Batteries AA / Mignon
- 4 rechargeable batteries (rechargeable with additional charger; REF 919 221)
- Battery pack, (REF 919 201) rechargeable via USB interface (starting with SN: NPF121600)

2.5. Operating environment

The following conditions are necessary to provide longevity and flawless functioning of the instrument.

- Operating temperature 0–50 °C
- Relative humidity: 90 %
- To ensure safe results, the instrument should not be stored in extremely dusty, damp or wet conditions (the instrument optics are IP67 proof and don't experience any sustained damage from such conditions).
- Keep the instrument surface, cuvette slot and any accessories clean and dry. Splashes or spilled substances should be wiped away and cleaned immediately (*see 8. Cleaning the photometer, page 30*)

Important note: Protect the instrument from extreme temperatures (e.g. radiators, direct sun light and other heat sources)

3. Usage

Due to unauthorized changes and modifications to the instrument without prior permission of authorized bodies, it is possible to revoke the users' operating permission.

As part of GLP, please familiarize yourself with all reagents used in your application. Read all product labels and MSDS before working with chemicals

To ensure proper functioning of all inbuilt safety measures, the instrument may, UNDER NO CIRCUMSTANCES, be used or installed in a different way than stated in the manual.

4. Setup

4.1. Unpacking the instrument

Take the instrument out of the box and check for damage that may have been caused during shipment. Please check if any parts listed in the package content are missing. If parts are missing or damaged, contact MACHEREY-NAGEL or your responsible distributor.

4.2. Operation with batteries

The PF-12 photometer runs on 4 AA alkali-manganese batteries (included in delivery). Rechargeable batteries, e.g. 4 Compit NiMH AA 2500 mAh can be used in the instrument and are available as accessories (REF 919 221). The charger and rechargeable batteries need to be ordered separately.

When changing batteries, we recommend to always exchanging the entire set. The battery compartment is at the bottom of the instrument. To change the batteries,

1. Turn the instrument around and place it on a soft surface.
2. Unscrew the screws of the battery compartment and take of the lid
3. Insert the 4 batteries in the respective slots. Make sure polarity of batteries and slots match.
4. Put the lid back on, tighten the screws and turn the instrument around.

Only remove the battery compartment when changing batteries!

4.3. Operation with battery pack

The PF-12 can be equipped with an internal battery pack, which can be charged via the USB port. You can order such a power pack separately and insert it into the instrument (as of SN: NPF121600). Detailed installation instructions come with the battery pack.

5. Operation

The photometer PF-12 has a 20-part keyboard:
10 number-keys and 10 keys with special functions

5.1. Keyboard

Key :

Keys  ... :

Keys  /  / :

Key :

Key :

Key :

Key :

Key :

Key :

Key :

Key :

Key :

On / Off switch

For entering digits, letters, special characters

Double-function keys, for sample number, sampling location (alphanumeric) and dilution

For entering the decimal point, to select special methods menu

Abort-key, return to method selection

Scroll-mode up key

Scroll-mode down key

Enter key, call up last result

Memory-key, call-up memory functions, store / not-store results

Extinction program-key, recall the extinction after measurement, display additional information

Manual zero adjustment

Measurement key, consecutive measurements are numbered continuously, Call up last test

5.2. Configuration

01.02.12 12:15

CONFIGURATION

→ (GB) ---> (?)
date / time
calibrate
memory

To get to the configuration mode enter    in the method selection mode (use  /  to select menu items). Press  to confirm and follow the instructions in the display.

1. Language

Once desired language appears, hold down  until the instrument switches back to configuration mode.

2. Date/Time

Navigate to date / time using  /  and adjust using the keyboard.

3. Calibration

Use  to turn on calibration. First the photometer automatically checks all filters with the lamp. For the second calibration step the photometer PF-12 requires a calibration tube (a clean *NANOCOLOR*® test tube filled with water). Confirm with .

4. Memory

Use / to turn the memory on or off and confirm by pressing .

5. Reaction time

Use / to turn the reaction time on or off and confirm by pressing .

6. Signal tone

Use / to turn the signal tone on or off and confirm by pressing .

7. Heading

Use the arrows / to select the desired text passage and add the heading with the alphanumeric keyboard.

8. Filter 8

Use to select the position of filter 8 and enter the wavelength of the additional special filter via the keyboard.

9. AUTO-OFF

Using / you can determine the photometer's auto-off time after the last measurement. Available times range from 5 to 20 min in 5-minute intervals(ex factory= 5 min), or you can disable the auto-off completely. Confirm your setting with .

10. SYSTEM-RESET

After an additional security confirmation the photometer performs a SYSTEM-RESET, thus restoring the state of delivery. The display then asks for a calibration. Confirm with key and proceed according to point 3.

11. PROG. UPDATE

In order to perform a program update, you must first install the update software and the new program file (www.mn-net.com) on your computer. Before you start the update procedure the photometer has to be connected to the computer via the USB-interface. Start the update program on the PC and set the photometer in the update mode using the scroll and confirmation key . After a successful update the photometer restarts.

5.3. Performance of **VISOCOLOR® ECO** tests

Instead of the **VISOCOLOR® ECO** measuring tube, a clean, dry 16 mm OD round tube is used. Zero is adjusted with a second round tube (depending on the test either an untreated sample or a reagent blind value).

01.02.12 12:15

Method: ---

Switch on photometer

The display „Method: _ _ _“ appears.

You can call up a **VISOCOLOR® ECO** Test test via the test number (5 - xx)

or

01.02.12 12:15

Method:

Basic functions
Tube tests
→ VISOCOLOR-ECO
Bio tests

with the keys / via a list.

01.02.12 12:15
5081 AMMONIUM 3
measure blank →N

Perform a zero measurement

Put the clean round tube with the blank value into the photometer and press . Remove round tube with blank value.

01.02.12 12:15
5081 AMMONIUM 3
measure sample →M

Perform a measurement

Put clean round tube with sample solution into the photometer, press  and read the result.

For the detailed procedure of the individual tests please follow the corresponding instruction leaflet. You can find detailed descriptions and pictograms on our homepage, www.mn-net.com, and on the software DVD.

5.4. Performance of **NANOCOLOR®** tube test

The characteristics of all our **NANOCOLOR®** tube tests are prestored in the PF-12.

01.02.12 12:15
Method: ---

Switch on photometer

The display “Method: _ _ _” appears.

You can call up a **NANOCOLOR®** tube test via the a test number (0 - xx)

or

01.02.12 12:15
Method:
Basic functions
→ Tube tests
VISOCOLOR-ECO
Bio tests

with the keys  /  via a list.

01.02.12 12:15
0091 LEAD 5
measure blank →N

Perform a zero measurement (for a few tests)

Put the clean round tube with the blank value into the photometer and press . Remove round tube with blank value

01.02.12 12:15
0032 AMMONIUM 3
measure sample →M

Perform a measurement

Put clean round tube with sample solution into the photometer, press key  and read the result.

For the detailed procedure of the individual tests please follow the corresponding instruction leaflet. You can find detailed descriptions and pictograms on our homepage, www.mn-net.com, and on the software DVD.

5.5. Basic photometric functions

01.02.12 12:15

Method:

- 901 FACTOR
- 902 STANDARD
- 903 EXTINCTION
- 905 TRANSMISSION

By calling up method 901 - 903 and 905 you can access basic photometric functions or select the required option in the scroll mode by using keys  /  and confirm with .

Method 901: Measurement with factor

1. Select the required wavelength using  /  and confirm the selection by pressing .
2. Enter factor via keyboard. Confirm with .
3. Perform measurement following the instructions in the display.
4. Press  to return to method selection. The factor will be erased.

Method 902: Measurement with standard

1. Select the required wavelength using  /  and confirm the selection by pressing .
2. Enter standard concentration via the keyboard. Confirm with .
3. Perform measurement following the instructions in the display.
4. Press  to return to method selection. The standard concentration will be erased.

Method 903: Measurement of extinction

1. Select the required wavelength using  /  and confirm the selection by pressing .
2. Perform measurement following the instructions in the display.
3. Press  to return to method selection.


Method 905: Measurement of transmission

1. Select the required wavelength using  /  and confirm the selection by pressing .
2. Perform measurement following the instructions in the display.
3. Press  to return to method selection.

6. Memory


In the standard mode of operation the memory is activated, indicated by a ✓ in the upper start display. The configuration menu allows to activate / deactivate memory completely.

For memory administration press  in the method selection mode. The display will read either,

01.02.12 12:15 ✓
MEMORY: 000

or, if any measuring values are stored, e.g.

01.02.12 12:15 ✓
MEMORY: 046/046
output ?

The memory administration offers several options for handling data. Using  /  you can select one of the four options.

Output

When you press  at this point, all stored data are transferred via the USB-Interface to a terminal program.

Display

When you press  at this point, you can display all data sets in sequence using . Using  allows you to return to the last-displayed value. By pressing  and confirming with  you can delete any data set.

Delete

When you press  at this point, and confirm deletion by pressing  once more, all data or the selected data are deleted from the memory.

Select

When you press  at this point, you access the selection mode.


The following options for selection are available:

Method / Location / Date / Time

7. Special functions


7.1. Special methods

The PF-12 allows you to define up to 20 special methods. The methods can be defined as linear (factor) or non-linear methods (up to the 4th degree). For defined methods, the same options as in the main program, such as sample location, dilution, memory, etc. are available.

Selection

Select special method from the main menu with .

Create new

Use  and e.g. 16 to create the new special method P16.


Now follows the menu for entering the method characteristics. To enter text, filter query, factor, measurement direction, measurement range, unit and dimension addition, you can enter letters and numbers via the key pad.

At the query <factor linear?>, you have two options. You can press  to enter a linear factor ($\beta = E \times F$). In this case, the decimal place of the factor also indicates the displayed decimal places.

Alternatively, you can press  /  to enter a polynomial up to the 4th degree ($\beta = F4 \times E4 + F3 \times E3 + F2 \times E2 + F1 \times E + F0$). Pressing button  before entering a polynomial factor results in a change of the algebraic sign. Once you have entered the polynomial factors, you can enter the number of decimal places in the display and the measurement direction (+ = increasing, - decreasing, +- both directions).

You can delete and re-enter an erroneous entry at any time by pressing . You can also abort programming a special method by pressing .

To run a special method, press  and enter the 2 digit number. The scroll buttons  /  offer the following options:

Execute

To select a previously programmed method, enter the 2 digit special method number, e.g. P18.

If the entered number does not refer to a special method, you can search for previously programmed methods by pressing  and then  / .

Press  to run the method.

01.02.12 12:15
P01 MACHEREY-NAGEL
edit ?

Edit

Press  /  to get to the edit sub-menu. Confirmation with  allows you to edit all previously made entries.

01.02.12 12:15
P01 MACHEREY-NAGEL
output ?

Export

Press  /  to get to the export sub-menu. Confirm with  to send all previous entries to a computer via the USB port.

01.02.12 12:15
P01 MACHEREY-NAGEL
delete ?

Delete

Press  /  to get to the delete sub-menu. Confirm with  and also confirm the following security question to delete the respective special method.

7.2. Sample numbering

To differentiate between samples, the PF-12 assigns a sample-number to every measurement (starting with 0001). Subsequent measurements are numbered accordingly. You have the option to change the sample-number after every measurement and to assign it to the current result. To do so, press  before the next measurement. An entry option appears to enter the 4 digit sample number via the key-pad.

Once you have entered the sample-number, it is displayed together with all other result data at the bottom of the display. The following measurements are then automatically numbered. Generally, the sample-number starts with 0001 whenever you select a method .

7.3. Identification of sample location

To differentiate between sampling locations, you have the option to enter a sample location after performing a measurement (up to 12 characters). Press button  to open the entry option for the sample location. Once you have entered the sample location, it is displayed together with all other result data at the bottom of the display. If you have entered additional sample locations before, you can press  /  to choose or edit the sample location from a list of up to 20 previously defined sample locations.

7.4. Sample dilution

During chemical analyses often samples have to be diluted to fit within a certain measuring range. To avoid manual calculations, you can enter the dilution into the PF-12, which will then re-calculate the result, taking the dilution into account.

Example:

Expected result: between 80 and 200 mg/L

Measuring range photometer: 0.1–10.0 mg/L

Necessary dilution (200 → 10 mg/L): at least 1:20

Recommended dilution: 1:25 to 1:50 to get within the middle of the measuring range

Dilutions are entered as 1 part sample plus x parts distilled water as to create distinct readings. You can enter dilutions between 1+1 and 1+999.

Press **[F3]** before the next measurement to open the entry option for the dilution. Once you have entered the dilution, it is displayed together with all other result data at the bottom of the display. The result is immediately re-calculated using the entered dilution.

Note: The entered dilution only impacts the current result. In the case of large dilutions, e.g. 1+999, it is possible that the result is displayed in a different dimension as to enable proper display in the screen.

8. Cleaning the photometer

You can clean the instrument's housing with a damp cloth. Splashes need to be wiped away immediately. You can clean and dry the cuvette slot using a cotton cloth.

Photometer and cuvettes need to be kept clean at any time. Cuvettes can be cleaned using a lint-free cloth.

9. Trouble shooting

Trouble shooting for problems arising in the field is limited to changing the batteries or the battery pack and calibrating the instrument.

The user may not perform any other maintenance activities. There are no other parts that may be exchanged or maintained. Opening the instrument housing terminates any warranty claim.

Disposal in accordance with EU directive 2002/96/EC

In compliance with EU directive 2002/96/EC MACHEREY NAGEL will accept back and dispose of the old instrument at no cost

ATTENTION!

It is not allowed to dispose of the instrument using municipal waste disposal services. Please talk to your local MACHEREY-NAGEL contact.

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